

AGE Reader



Non-invasive assessment of cardiovascular risk

Introduction



- Cardiovascular risk
- Measuring CV risk
- AGE Reader
- Clinical validation

Cardiovascular risk

- 70% of diabetics will die of CV event
- 5 out of 7 drugs (type 2 diabetic)
- Events are costly
- Difficult to assess

Measuring CV risk

Now in the clinic:

- Risk scores (UKPDS, SCORE, Framingham)
 - Intima-media thickness
 - Coronary calcification (CT scan)
-
- The AGE Reader gives
 - quickly and non-invasively
 - reliable extra information to help the physician to improve therapy


AGEs



- Advanced Glycation Endproducts
 - Glycated proteins
- AGE values rise in everyone when getting older.
 - But, much quicker in diabetes
- AGEs play a central role in the development of diabetes and cardiovascular risk.


AGEs

- AGE formation





AGEs

- AGEs causes major damage to tissue and cells
 - Irreversible malfunctioning proteins
 - Extreme long $\frac{1}{2}$ life (many years)
 - In particular in diabetic patients



AGE Reader



AGE Reader



■ Measurement report




Measurement report

Number: 01
Name: Raymon
Gender: Male
Age: 24

Measurement Results

AF 1.3

Measurement setting: Triple Measurement
Measured on: 2-6-2010 10:48



AGEs

DiagnOptics AGE Reader is a medical device to estimate cardiovascular risk. The AGE Reader non-invasively assesses the accumulation of advanced glycation endproducts (AGEs) in the skin using fluorescence of ultraviolet light. AGEs play a pivotal role in the development of chronic complications of diabetes and other common conditions. The amount of AGEs in tissues serves as an important risk predictor of such complications.

Breakthrough publication



Diabetologia (2009) 52:789–797
DOI 10.1007/s00125-009-1308-9

ARTICLE

Skin autofluorescence provides additional information to the UK Prospective Diabetes Study (UKPDS) risk score for the estimation of cardiovascular prognosis in type 2 diabetes mellitus

H. L. Lutgers · E. G. Gerrits · R. Graaff · T. P. Links · W. J. Shuiter · R. O. Gans · H. J. Bilo · A. J. Smit

Received: 24 October 2008 / Accepted: 30 January 2009 / Published online: 10 March 2009
© The Author(s) 2009. This article is published with open access at Springerlink.com

Abstract

Aims/hypothesis The UK Prospective Diabetes Study (UKPDS) risk engine has become a standard for cardiovascular risk assessment in type 2 diabetes mellitus. Skin autofluorescence was recently introduced as an alternative tool for cardiovascular risk assessment in diabetes. We investigated the prognostic value of skin autofluorescence for cardiovascular events in combination with the UKPDS risk engine in a cohort of patients with type 2 diabetes managed in primary care.

H. L. Lutgers (✉) · R. O. Gans · H. J. Bilo · A. J. Smit
Department of Internal Medicine,
University Medical Center Groningen,
Hanzeplein 1,
9700 RB Groningen, the Netherlands
e-mail: h.l.lutgers@int.umcg.nl

E. G. Gerrits · H. J. Bilo
Diabetes Center, Isala Clinics,
Zwolle, the Netherlands

R. Graaff
Department of Biomedical Engineering,
University Medical Center Groningen,
Groningen, the Netherlands

T. P. Links · W. J. Shuiter
Department of Endocrinology,
University Medical Center Groningen,
Groningen, the Netherlands

T. P. Links · R. O. Gans · H. J. Bilo · A. J. Smit
Department of Medicine, University of Groningen,
Groningen, the Netherlands

Methods Clinical, UKPDS risk engine and skin autofluorescence data were obtained at baseline in 2001–2002 in the type 2 diabetes group ($n=973$). Follow-up data concerning fatal and non-fatal cardiovascular events (primary endpoint) were obtained till 2005. Patients were classified as ‘low risk’ when their 10 year UKPDS risk score for fatal cardiovascular events was <10%, and ‘high risk’ if >10%. Skin autofluorescence was measured non-invasively with an autofluorescence reader. Skin autofluorescence was classified by the median (i.e. low risk < median, high risk > median).

Results The incidence of cardiovascular events was 119 (44 fatal, 75 non-fatal). In multivariate analysis, skin autofluorescence, age, sex and diabetes duration were predictors for the primary endpoint. Addition of skin autofluorescence information to that from the UKPDS risk engine resulted in re-classification of 55 of 203 patients from the low-risk to the high-risk group. The 10 year cardiovascular event rate was higher in patients with a UKPDS score >10% when skin autofluorescence was above the median (55.8% vs 38.9%).


Conclusions/interpretation Skin autofluorescence provides additional information to the UKPDS risk engine which can result in risk re-classification of a substantial number of patients. It furthermore identifies patients who have a particularly high risk for developing cardiovascular events.

Keywords Advanced glycation end-products · Complications · Fluorescence · Type 2 diabetes · UKPDS risk score



- Two important conclusions:
 - AGE measurement is best single predictor of (CV) mortality (after age)
 - When compared with whole risk engine:
 - independent
 - 27% re-classification

Additive predictive value skin AF on top of UKPDS risk score in type 2 diabetes (n=967)



News Release

FOR IMMEDIATE RELEASE

Contact: Rachel Morgan, ADA
(703) 549-1500 ext. 2290

New Tool Can Help Predict Diabetes Complications

*Other Studies This Month: Coffee Helps Prevent Diabetes;
When You Get Diabetes (Middle v. Old Age) Matters*

(Alexandria, VA) – A noninvasive tool that measures the skin’s autofluorescence could help doctors determine whether people with diabetes are beginning to develop serious complications, according to a study published in the November issue of *Diabetes Care*.

Researchers in the Netherlands found that illuminating a patient’s lower arm with a fluorescent tube accurately reflects vascular damage caused by the accumulation of advanced glycosylation end products (AGEs). AGEs are produced in the body when glucose links with protein. They play a role in damaging blood vessels, which can lead to complications, such as nerve damage.

Previous studies have shown that AGEs have fluorescent properties. This study confirmed that those properties could be measured by illuminating the skin, and that high levels of autofluorescence were associated with more severe diabetes complications, such as neuropathy, retinopathy and cardiovascular problems.

“With this tool, doctors could easily check people with diabetes in an outpatient clinic setting to see whether they may already be developing dangerous complications,” said lead researcher Dr. Helen Lutgers, of the University Medical Center in Groningen, the Netherlands. “The sooner complications are detected, the better the chance of preventing progression of damage.”

The technology used in this study is currently commercially available in Europe. Until FDA approval is obtained, its availability in the USA is restricted to experimental use only.

To reach Dr. Lutgers, email H.L.Lutgers@isala.nl or phone 31-5036-10096.

Case Report 1

AGE Reader and CV risk
in clinical practice

Case Report

UKPDS Risk Engine v2.0

Input





Age now: years HbA1c: %
Diabetes duration: years Systolic BP: mm Hg
Sex: Male Female Total cholesterol: mmol/l
Atrial fibrillation: No Yes HDL cholesterol: mmol/l
Ethnicity:
Smoking:

Number of values*

HbA1c:
Systolic BP:
Total cholesterol:
** used to adjust for regression dilution*
Units: mmol/l mg/dl

Output

10 year risk 0 15 30 100

| | | |
|---------------|------------------------------------|--|
| CHD: | <input type="text" value="15.9%"/> |  |
| Fatal CHD: | <input type="text" value="8.7%"/> |  |
| Stroke: | <input type="text" value="4.5%"/> |  |
| Fatal stroke: | <input type="text" value="0.7%"/> |  |

Adjusted for regression dilution

Risk interval

Risk over next years

Ideal doctor:
“when do I
need to
intensify
treatment?”

→ Risk
often in
‘grey’
zone

Ideal doctor knows: guidelines tell me to weigh other risk factors, but

- Metabolic syndrome, CRP, homocysteine: only marginal additive (to UKPDS RE) improvement in prediction of CV events
- Microralbuminuria, renal dysfunction, hsCRP: require repeated sampling because of variability, medication, intermittent disease etc.
- IMT or coronary calcification score not generally available, expensive, radiation exposure. Still > 4-fold increase coronary CT in US over last decade (NEJM 2009)

Consequences of current guidelines in daily clinical practice

- Clinicians aware of weaknesses of risk scores
- Use UKPDS RS requires all 10 items and risk score result to be available at moment of patient visit to make decisions. In practice not so easy ...
- Therefore, risk engines not used at all in majority!

- Similar situation for use of additional RF
- Then: to be on the safe side, doctor decides to intensify treatment

Case Report 2

CV risk in clinical practice

+

AGE Reader measurement

Clinical practice

- AGE Reader in clinical practice

UKPDS Risk Engine v2.0

Input

Age now : years HbA1c : %
 Diabetes duration : years Systolic BP : mm Hg
 Sex : Male Female Total cholesterol : mmol/l
 Atrial fibrillation : No Yes HDL cholesterol : mmol/l
 Ethnicity :
 Smoking :

Output

| | 10 year risk | 0 | 15 | 30 | 100 |
|----------------|------------------------------------|---|----|----|-----|
| CHD : | <input type="text" value="11.0%"/> | | | | |
| Fatal CHD : | <input type="text" value="7.6%"/> | | | | |
| Stroke : | <input type="text" value="5.9%"/> | | | | |
| Fatal stroke : | <input type="text" value="1.0%"/> | | | | |

Adjusted for regression dilution

AF value = 3.6.


Additive info with AGE reader:

If AF = 3.6 → median life expectancy this patient $25/2.62 = 9.5$ yr

If the same patient has an **AF = 2.4**, → Median life expectancy 14.4 yr


Extension to T1DM

- Chabroux et al. (2010):
- French study in Lyon region in T1DM
- Confirms associations between AGE reader levels and nephropathy and neuropathy complications



Added proof in T2DM;

- AURORA study:
- Cross-sectional study in 566 T2DM patients controlled in 5 Dutch hospitals: sicker, higher HbA1c, more comorbidity
- Relation SAF with complications similar to previous study in primary care, better controlled patients
- SAF adds again to UKPDS!



Independent contribution SAF to UKPDS RS for existing CVD

Renal failure

- CV complications and mortality in renal failure:
- In dialysis patients: 2 new, Japanese and Dutch , cohorts. Now studies available in 3 different cohorts on predictive value in hemodialysis patients
- More importantly: association SAF with CV complications extends also to CKD 3 (GFR 30-60 ml/min)(McIntyre, RRID; CJASN 2011):

**Therapeutic implications more important
in this group !**

Skin Autofluorescence and the Association with Renal and Cardiovascular Risk Factors in Chronic Kidney Disease Stage 3

Natasha J. McIntyre, Richard J. Fluck,* Christopher W. McIntyre,*† and Maarten W. Taal**

Clin J Am Soc Nephrol 6: ●●●–●●●, 2011. doi: 10.2215/CJN.02420311

Conclusion

We have shown that tissue AGE, measured as SAF, is independently associated with multiple traditional and nontraditional risk factors for CKD progression and CVE in community-based participants with CKD stage 3. SAF measurement may therefore represent a clinically useful, noninvasive method for assessing renal and cardiovascular risk in participants with CKD. Planned follow-up for up to 10 years will evaluate this hypothesis in the RRID cohort.

USPs AGE Reader

- Clinically proven
- Non-invasive
- Reproducible
- No consumables
- Quick
- Easy to perform

Thank you for your
attention

