




Name: John Doe
Year of birth: 1960
Gender: Male
Measurement date: 18-03-11

Result:

Based on the measurement result and the data below you are at increased risk for having (pre-)diabetes. We advise you to visit your doctor for an official diabetes diagnostic test.



*The normal range is normalized and corrected for calendar age.

Auto-Immune disease: Negative
Hospitalized (<6 months): Negative
Kidney disease: Negative

Diab-spot is a medical device, developed by DiagnOptics, to estimate the risk for (pre-)diabetes. Diab-spot non-invasively assesses the accumulation of Advanced Glycation Endproducts (AGEs) in the skin using fluorescence of ultraviolet light. The amount of AGEs in tissues serves as an important risk predictor of (pre-)diabetes. Diab-spot combines the AGE measurement with some clinical characteristics to calculate the Diab-spot test result.